Poster 950



Evaluation of the Performance of the Target Enriched Multiplex Polymerase Chain Reaction (TEM-PCRTM) Gastrointestinal Panel vs the Cepheid Xpert® C. difficile/Epi Assay for the **Detection of Toxigenic** *Clostridium difficile*

Khushdeep Chahal, MD1, Esmeralda Gutierrez-Asis, MD1, Farrah Ibrahim, MD FACP1, Elena Grigorenko, PhD2, Stefan Brzezinski2, Donald Stalons, PhD, D(ABMM), MPH2, Ali Hassoun, MD FIDSA FACP3

Internal Medicine, University of Alabama at Birmingham Huntsville Regional Medical Campus, Huntsville, AL, 2Diatherix Laboratories, LLC, Huntsville, AL, 3University of Alabama School of Medicine-Huntsville Campus, Huntsville, AL

Khushdeep Chahal, MD chahalk@uab.edu Ph. No: 256-551-4557 Abstract No: 53254 Stefan Brzezinski stefan.brzezinski@diatherix.com

Abstract:

Background: Clostridium difficile (CD) is one of the most common healthcare associated infections (HAIs). While the gold standard for diagnosis of CD is cytotoxicity assay and toxigenic culture, there are a number of molecular assays for rapid CD detection. This study evaluated the performance of an in-house TEM-PCR™ GI panel with the Cepheid Xpert® C. difficile/Epi assay using clinical samples collected in hospital setting.

Methods: Stool samples of patients with suspected CD infections were tested by TEM-PCR GI panel and with Cepheid Xpert® assay. CD detection is the part of the Diatherix multiplex GI panel which can simultaneously detect 3 viral, 8 bacterial, 2 protozoan, and 2 toxins in a single specimen. The Cepheid assay detects the toxin B gene and the O27/NAP1/B1 strain. Two TagMan® assays targeting the toxin B and to gene were developed as an alternative approach for confirmation of results. Demographics data, risk factors, and other clinically relevant information was

Results: We enrolled 80 patients with mean age of 62 years (TaqMan® assays done in 80). 48% females. Known risk factors included: rec Results: We enrolled 80 patients with mean age of 62 years (TaqMan09 assays done in 80). 48% temales. Known risk ractors includer recurrent CD infection (15%), recent hospitalization (55%), proton pump inhibitor use (45%), gastrointestal disorders (39%) and immunosuppression (19%). There was 98.8% concordance between TEM-PCR CD toxin B, TaqMan09 typene, and Cepheid Xpert®. Cepheid NAPI/B1 was positive in 8% of patients. EIN-PCR detected positive binary toxin in 14% of patients. All NAPI/B1 positive specimers were positive for binary toxin with TEM-PCR. 17% of patients with positive CD by TEM-PCR had co-detection of other targets (25% Norovirus, 50% Rotavirus and 25% EPCE). Those patients with co-detection, 75% were female and 75% more than 55 years old. Othool studies done in 63% of patients, 44% had occult blood with 71% positive, 39% had cultures with 6% positive for Gram positives, and 40% had O&P (All were negatve).

Conclusion: The performance of TEM-PCR GI panel was comparable to Cepheid assay for CD detection. Further studies are needed to asses the clinical utility of TEM-PCR GI panel to detect multiple pathogens and how this information can be used to improve patient care

Note: abstract updated to reflect the inclusion of additional patients for negative predictive value calculations

Introduction:

- · C. difficile infection (CDI) is a leading cause of nosocomial infectious diarrhea and is emerging as a community acquired
- pathogen in groups that were previously considered low risk.⁵
 Symptoms of CDI can be caused by other GI pathogens that may be missed if only a *C. difficile* assay is used.
- This study sought to ask the question if a multiplex PCR panel like the Diatherix TEM-PCR™ GI Panel could detect CDI as well as the FDA-cleared Cepheid Xpert® C. diff/Epi assay while giving the added benefit of multiplex pathogen detection.
- Patients in this study were inpatients at a local hospital with symptoms suspected to be associated with Clostridium difficile infection.
- Stool specimens were split and tested with the Cepheid C. diff/ Epi FDA cleared assay and by Diatherix Laboratories on their TEM-PCR™ Gastrointestinal Panel.
- All samples were also tested with a C. difficile binary toxin end-point assay, a C. difficile toxin B TaqMan® assay, and a $\textit{C. difficile} \ \ \text{TaqMan} \\ \text{\textcircled{\textbf{B}} assay that identifies generalized } \textit{C. difficile} \ \ \text{regardless of toxin presence. A combination of three}$ confirmatory assays were used and a sample was considered positive for *C. difficile* if two out of three of these assays were
- Samples were blinded when tested with the Diatherix panel

Materials & Methods:

Specimens were extracted using a modified Omega Bio-Tek, Inc. Mag-Bind® kit (Norcross, GA) with additional reagents from MO BIO Laboratories, Inc. (Carlsbad, CA) on a ThermoFisher Scientific, Inc. KingFisher™ Flex (Waltham, MA) semi -automated platform.

Extracted DNA and RNA samples were amplified using the Target Enriched Multiplex PCR (TEM-PCR*) Gastrointestinal Panel.TEM -PCR** is a highly multiplexed, nested, end-point PCR technique covered by US patent 7,851,148 B2. The Gastrointestinal Panel contains 59 target specific primers which can identify 13 different gastrointenstinal pathogens. The end-point PCR products generated by TEM-PCR** were detected on a custom microarray (Microarrays, Inc., Huntsville, AL).

PCR and Detection for TagMan® gPCR Confirmation

R after Detrection for Lagination GPC Accommission of April 2007 (April 2007) and Apri cycles in a BioRad CEX96 (BioRad, Hercules, CA), PCR was done in triplicate.




Figure 1. TEM-PCR™ Scheme.

Figure 1. IEM-PCR Scheme.
Low concentrations of nested gene-specific primers (Fo-forward out; Fi-forward in; Ri-reverse in; and Ro-reverse out) are designed to enrich the genetic targets during the initial PCR cycles. Later in the procedure, a pair of universal SuperPrimers (FS and RS) are used to amplify all targets. The reverse SuperPrimer is labeled with bloin. Target-specific biotinylated PCR products are detected with a complimentary detection serves which is convenient countries. detection probe which is covalently coupled to a glass

Diatherix TEM-PCR™ Gastrointestinal Panel

Campylobacter jejuni Giardia lamblia Vibrio parahaemolyticus Enteropathogenic E. coli (EPEC)

Enterotoxigenic E. coli (ETEC) Enterohemorrhagic E. coli (EHEC) Shiga-toxin 1 (stx1) Shiga-toxin 2 (stx2) Shigella/ Enteroinvasive E. coli (EIEC) Salmonella enterica Adenovirus 40/41

Results:

Table 1. Oligonucleotide sequences used for confirmatory testing of toxigenic C. difficil

	Target Gene	Oligonucleotide	Sequence (5'-3')	Position	Amplicon size (bp)	GenBank Accession No.	
tcdB	C. difficile	Forward	TGATTGCAGTTGTAGCTGTTGTTAAA	3621-3646			
	toxin B	Reverse	CGAGTGACCCATTATTAAGACAAGAA	3563-3588	58	AF217292	
		Probe	FAM-TTACTGCCATTATACCTATCTT-MGB	3598-3619			
tpi	C. difficile	Forward	AAGCATTAGAAGTAGGAATAGACCCAAT	315-342			
	triosephosphate	Reverse	TTTAGTTTTTCCAGCTTCTCTTTGTTC	365-391	76	AY700149	
	isomerase	Probe	FAM-TTATGTGTTGGAGAAACT-MGB	344-361			

 Table 2. Comparison of detection results of C. difficile between Cepheid, TEM-PCR, and confirmatory methods, such
 as a combination of TagMan® toxB and toi assays and TEM-PCR based detection of binary toxin

Assay	Cepheid C. diff / Epi		Assay performance (95% confidence interval)					
Assay	Detected	Not Detected	Sensitivity (%)	Specificity (%)	PPV (%)	NPV (%)		
Diatherix TEM-PCR C. diff toxin B								
Detected	24	0						
Not Detected	1	55	96 (77.7-99.8)	100 (91.9-1)	100 (82.8-1)	98.2 (89.2-99.9)		
Total	25	55						
Confirmatory m	ethods							
Detected	23	0						
Not Detected	2	55	92 (72.4-98.6)	100 (91.9-1)	100 (82.2-1)	96.5 (86.8-99.4)		
Total	25	55						

CI calculated with efficient score method with continuity correction

Table 3. Results of the confirmatory assays of specimens positive for toxigenic *C. difficile* by Cepheid results including an end-point binary toxin assay (*cdtA*, *cdtB*), a toxin B TaqMan® qPCR assay, and a *tpi* TaqMan® qPCR assay for general *C. difficile*. The qPCR assays are less sensitive than the TEM-PCR assay and therefore a combination of three


Cepheid C. difficile Result	TEM-PCR Result	Binary toxin?	ToxB Taqman Result	TPI Taqman Result	Concordance?	Decision Criteria
DETECTED	DET	YES	DET	DET	YES	BT/TOB/TPI
DETECTED	DET	NOT DET	DET	DET	YES	TOB/TPI
DETECTED	DET	NOT DET	NOT DET	DET	NO	TPI
DETECTED	DET	YES	DET	DET	YES	BT/TOB/TPI
DETECTED	DET	YES	DET	DET	YES	BT/TOB/TPI
DETECTED	DET	YES	NOT DET	DET	YES	BT/TPI
DETECTED	DET	NOT DET	DET	DET	YES	TOB/TPI
DETECTED	DET	NOT DET	DET	DET	YES	TOB/TPI
DETECTED	DET	NOT DET	DET	DET	YES	TOB/TPI
DETECTED	DET	YES	DET	DET	YES	TOB/TPI
DETECTED	DET	YES	DET	DET	YES	BT/TOB/TPI
DETECTED	DET	NOT DET	DET	DET	YES	BT/TOB/TPI
DETECTED	DET	NOT DET	DET	DET	YES	TOB/TPI
DETECTED	DET	NOT DET	DET	DET	YES	TOB/TPI
DETECTED	DET	YES	NOT DET	DET	YES	BT/TPI
DETECTED	DET	YES	DET	DET	YES	BT/TOB/TPI
DETECTED	DET	NOT DET	DET	DET	YES	TOB/TPI
DETECTED	DET	NOT DET	DET	DET	YES	TOB/TPI
DETECTED	DET	NOT DET	DET	DET	YES	TOB/TPI
DETECTED	DET	YES	DET	DET	YES	BT/TOB/TPI
DETECTED	DET	YES	DET	DET	YES	BT/TOB/TPI
DETECTED	DET	NOT DET	DET	DET	YES	TOB/TPI
DETECTED	DET	YES	DET	DET	YES	BT/TOB/TPI
DETECTED	DET	NOT DET	DET	DET	YES	TOB/TPI
DETECTED	NOT DET	NOT DET	NOT DET	MOT DET	NO	N/A

Table 4. TaqMan® assay specificity using high titers of organisms tested in duplicate

Sample	Conc.	tcdB Cq		
NTC		ND	ND	Gardner
Aeromonas hydrophila 35654	1.5e8 cfu/mL	ND	ND	Haemop
Atopobium vaginae BAA-55	1e7 cfu/mL	ND	ND	Hafnia :
Bacillus cereus 10876	1.5e8 cfu/mL	ND	ND	Helicob
Bacteroides thetaiotaomicron 29741	1.5e8 cfu/mL	ND	ND	Klebsie
Bifidobacterium longum/E. rectale 35183	1.5e8 cfu/mL	ND	ND	Lactoba
Campylobacter coli 43133	1.5e8 cfu/mL	ND	ND	Listeria
Campylobacter jejuni BAA-1153	1.5e8 cfu/mL	ND	ND	Morgan
Candida albicans 11006	1.5e8 cfu/mL	ND	ND	Peptos
Candida glabrata 32554	1.5e8 cfu/mL	ND	ND	Pleison
Chlamydia trachomatis 0801775	1e7 ifu/mL	ND	ND	Prevote
Citrobacter freundii 8090	1.5e8 cfu/mL	ND	ND	Proteus
Clostridium difficile 43596	1.5e8 cfu/mL	21.6 ^A	21 ^A	Pseudo
Clostridium difficile 700057 (non-toxigenic)	1.5e8 cfu/mL	ND	20.8 ^A	Salmor
Clostridium novyi 27606	3e8 cfu/mL	ND	ND	Salmor
Clostridium septicum 8065	3e8 cfu/mL	ND	ND	Serratia
Clostridium sordellii 9714	1.5e8 cfu/mL	22.3	22.3	Shigella
Enterobacter aerogenes 13048	1.5e7 cfu/mL	ND	ND	Shigell
Enterobacter cloacae 13047	1.5e7 cfu/mL	ND	ND	Staphy
Enterococcus faecalis 700802	1.5e7 cfu/mL	ND	ND	Staphy
Enterococcus faecium 700221	1.5e7 cfu/mL	ND	ND	Strepto
Eschericia coli O111:H8 BAA-2217 (EHEC)	1.5e8 cfu/mL	ND	ND	Strepto
Eschericia coli O157:H7 43895 (EHEC)	1.5e8 cfu/mL	ND	ND	Veillone
Eschericia coli O29:NM 43892 (EIEC)	1.5e8 cfu/mL	ND	ND	Vibrio v

The tcdB and tpi assays were found to cross-react with the closely related organism, Clostridium sordellii. C. sordellii produces toxins very similar to C. difficile toxin A and B. Antibodies to C. sordellii toxins will bind C. difficile toxins. C. sordellii is uncommon in humans. 1.2 A Expected result

Results (continued)



Conclusions:

- The C. difficile toxin B target in the Diatherix TEM-PCR™ Gastrointestinal Panel was determined to be 98.8% in agreement with the FDA cleared Cepheid C. diff/ Epi real-time PCR assay that targets the toxin B gene (tcdB, binary toxin gene, and hypervirulence associated tcdC deletion
- One sample was positive with the Cepheid test but was negative for TEM-PCR toxin B, TaqMan® toxin B, TaqMan® tpi, and TEM-PCR binary toxin suggesting a false positive Cepheid result or a mislabeling of the specimen.

 The performance of Diatherix Gastrointestinal Panel for CDI detection is comparable to the FDA cleared Cepheid C. diff/
- Epi test while offering simultaneous detection of 12 clinically relevant GI other pathogens.
- Norovirus, rotavirus, and EPEC were found to be the most common co-detections with C. difficile in the study population using large dataset of patient samples with suspected gastroenteritis submitted to Diatherix.
- As expected, recent antibiotic use was found to be the most prevalent risk factor for samples positive for CDI detection.

- Lyerly, D. M., Krivan, H. C., Wilkins, T. D. (1988). Clostridium difficile: Its Disease and Toxins. Clinical Microbiology Reviews. Jan. p. 1-18.
- Rovida, F., Campanini, G., Adzasehoun, K.M., Sarasini, A., Baldanti, F. Molecular Detection of Gastrointestinal Viral Infections in Hospitalized Patients. I Infectious Disease. 77 (3) 231-235
- Borrows, C.L., Turner, P.C. Seasonal Screening for Viral Gastroenteritis in Young Children and Elderly Hospitalized Patients: Is It Worthwhile? Journal of Hospital Infection. 87
- Ghose C. Clostridium difficile infection in the twenty-first century. Emerging Microbes & Infections. 2013;2(9):e62-. doi:10.1038/emi.2013.62

ents: Thanks to the Diatherix Clinical laboratory and everyone involved for processing the specimens. Thanks to Matthew Huff for collecting the prevalence data.

